

Related concepts

Linear expansion, volume expansion of liquids, thermal capacity, lattice potential, equilibrium spacing, Grüneisen equation.

Principle and task

The volume expansion of liquids and the linear expansion of various materials is determined as a function of temperature.

Equipment

Dilatometer with clock gauge	04233.00	1
Copper tube for 04231.01	04231.05	1
Aluminium tube for 04231.01	04231.06	1
Tube, quartz for 04231.01	04231.07	1
Immersion thermostat A100	46994.93	1
Cooling coil f.A100	46994.01	1
Accessory set for A100	46994.02	1
Bath for thermostat, Makrolon	08487.02	1
Lab thermometer, -10...+100C	38056.00	1
Rubber tubing, i.d. 6 mm	39282.00	2
Syringe 1ml, Luer, 10 pcs	02593.03	1
Cannula 0.6×60 mm, Luer, 20 pcs	02599.04	1
Measuring tube, l 300 mm, NS19/26	03024.00	2
Wash bottle, plastic, 250 ml	33930.00	1
Flask, flat bottom, 50 ml, IGJ19/26	35810.01	2

Glass beaker, tall,	100 ml	36002.00	1
Ethyl acetate	250 ml	30075.25	1
Glycerol	250 ml	30084.25	1
Olive oil, pure	100 ml	30177.10	1
Laboratory balance w. RS 232, 310 g		45025.93	1


Problems

1. To determine the volume expansion of ethyl acetate (C₄H₈O₂), methylated spirit, olive oil, glycerol and water as a function of temperature, using the pycnometer.
2. To determine the linear expansion of brass, iron, copper, aluminium, duran glass and quartz glass as a function of temperature using a dilatometer.
3. To investigate the relationship between change in length and overall length in the case of aluminium.

Set-up and procedure

1. The volume of the pycnometer is determined and the scale calibrated by weighing it empty and then filled with distilled water.
- The pycnometer, filled with the liquid to be measured, is brought to temperature in the water bath (thermostat). The change in volume is read from the scale on the tube built into its stopper.

Fig. 1: Experimental set-up for measuring thermal expansion.



2. The connecting tube to the thermostat is removed and the dilatometer is connected to the water circuit instead. Keep the feed and discharge lines as far away from the dilatometer as possible so that its body will not heat up. Clamp on the measuring tube, set the scale on the dial gauge to "0" and measure the expansion as a function of the temperature.

There is so little expansion in the case of duran glass and quartz glass that the heating and expansion of the dilatometer body as a result of radiation and conduction falsifies the measurement considerably. In this case, therefore, the measurement is started at the highest temperature (80°C) and the hot water in the bath replaced with cold tap water.

As the temperature changes very quickly with this method, the temperature of the dilatometer body remains constant. Only two values are measured.

3. In the case of aluminium, expansion is measured at three different rod lengths. The rod can be clamped in various places for this.

Theory and evaluation

An increase in temperature T causes the vibrational amplitude of the atoms in the crystal lattice of the solid to increase. The potential curve (Fig. 2) of the bonding forces corresponds only to a first approximation to the parabola of a harmonic oscillation (dotted line); generally it is flatter in the case of large interatomic distances than in the case of small ones. If the vibrational amplitude is large, the centre of oscillation thus moves to larger interatomic distances. The average spacing between the atoms increases, as well as the total volume V (at constant pressure p).

$$\alpha = \frac{1}{V} \cdot \left(\frac{\partial V}{\partial T}\right)_p \quad (1)$$

is called the volume expansion coefficient; if we consider one dimension only, we obtain the coefficient of linear expansion

$$\alpha_1 = \frac{1}{l} \cdot \left(\frac{\partial l}{\partial T}\right)_p \quad (2)$$

where l is the total length of the body.





Fig. 2: Potential curve as a function of the interatomic spacing r .

Fig. 3: Relationship between volume V and temperature ϑ of: a) ethyl acetate, b) methylated spirit, c) olive oil, d) glycerol and e) water.



1. A rise in the temperature causes a greater thermal agitation of the molecules in a liquid and therefore an increase in its volume (water between 0 and 4°C is an exception to this, however).

The coefficient of expansion of olive oil and water depends on temperature. Measured values at 20°C are:

	$\alpha/10^{-3}K^{-1}$
Water	0.20
Glycerol	0.50
Olive oil	0.72
Methylated spirit	1.11
Ethyl acetate	1.37

2. Fig 4 shows that the length increases approximately linearly with the temperature in the temperature range observed. Since the changes in length

$$\Delta l = l - l_0$$

are small compared with the original length l_0 , we can say

$$\alpha_1 = \frac{\Delta l}{l_0} \cdot \frac{1}{\Delta \vartheta} \quad (3)$$

and thus

$$l = l_0 [1 + \alpha_1(\vartheta - \vartheta_0)] \quad (4)$$

where ϑ_0 is the initial temperature.

The coefficients of linear expansion measured are:

	$\alpha_1/10^{-3}K^{-1}$
Aluminium	2.2
Brass	1.8
Copper	1.6
Steel	1.1
Duran glass	0.32
Quartz glass	0.046

The coefficient of expansion of steel and aluminium depends on the composition of the metal used.

3. If the temperature changes $\Delta\vartheta$ are not too large, the change in length Δl is proportional to the original length l_0 (See (3)).





Fig. 4: Relationship between length l and temperature ϑ , for a) aluminium, b) brass, c) copper, d) steel, e) duran glass, f) quartz glass ($l_0 = 600$ mm)

Fig. 5: Change in length Δl as a function of the original length l_0 for aluminium at $\Delta\vartheta = 15K$.



Note

The Grüneisen equation

$$\frac{\alpha}{C_p} = \gamma \cdot \frac{\kappa}{V} \tag{5}$$

where

$$\kappa = -\frac{1}{V} \left(\frac{\partial V}{\partial p} \right)_T$$

is the compressibility and

$$C_p = \left(\frac{\partial U}{\partial T} \right)_p$$

is the thermal capacity of the solid ($U =$ internal energy), signifies a relationship between the mechanical and thermal properties of a solid.

The Grüneisen parameter γ is defined by the change in the frequency ν of lattice vibration with volume:

$$\frac{\Delta \nu}{\nu} = -\gamma \frac{\Delta V}{V}$$

and can be calculated from macroscopic quantities in accordance with (5).